



Preparation of *trans*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride

Trans-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride is a dark green crystalline substance very well soluble in water and almost insoluble in common organic solvents. In aqueous solution it decomposes stepwise:




By evaporating a neutral solution at 90 – 100 °C the *trans*-isomer is changed to a racemic mixture of the Λ -*cis*-isomer and Δ -*cis*-isomer.

The *trans*-isomer may be prepared by oxidation of cobalt(II) chloride solution containing ethane-1,2-diamine. As an oxidizing agent gaseous oxygen may be used. After adding hydrochloric acid, heating and evaporating, a solid adduct *trans*-[CoCl₂(en)₂]Cl · HCl is formed.



By careful and moderate thermal decomposition of the adduct the *trans*-isomer is prepared.



Work

According to the following procedure prepare *trans*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride.

Chemicals

- cobalt(II) chloride hexahydrate, CoCl₂ · 6 H₂O, dark violet crystalline substance,
- ethane-1,2-diamine (ethylene diamine), NH₂CH₂CH₂NH₂, colourless oily liquid, $\rho = 0,898 \text{ g cm}^{-3}$,
- hydrochloric acid, HCl, concentrated, $w(\text{HCl}) = 0,36$,
- ethanol, CH₃CH₂OH, denatured spirit,
- diethyl ether, (C₂H₅)₂O.

Procedure

Dissolve 5,0 g of finely powdered cobalt(II) chloride hexahydrate in 20 cm³ of water. Add 2,8 cm³ (= 2,5 g) of ethane-1,2-diamine. Pour the final solution to a gas-washing bottle and connect it to a water pump (Fig. 1). Bubble the mixture with oxygen for about one hour. Put the reaction mixture to a beaker with hot water to speed up the oxidation reaction. Once the oxidation is finished, transfer the solution from the gas-washing bottle to an empty beaker. Under continuous stirring add 10 cm³ of concentrated hydrochloric. Bring the mixture in the beaker to a boil and evaporate *cca* half of its volume. Then, cool the mixture with ice. After cooling to room temperature, filter out the green adduct *trans*-[CoCl₂(en)₂]Cl · HCl on a Büchner funnel. Wash the adduct with ethanol and diethyl ether.

Put the adduct into an evaporating dish in thin layer, and heat it moderately over a ceramic wire gauze at temperature 110 – 125 °C for about 30 minutes. Eliminating gaseous hydrogen chloride a green crystalline *trans*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride will be prepared.





Fig. 1 Apparatus for preparation of *trans*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride.

1 – reactor heated by hot water,
2 – safety vessel.

Preparation of *cis*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride

Cis-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride is a violet crystalline substance very well soluble in water and almost insoluble in common organic solvents.

The procedure for preparation of the *cis*-isomer is based on stepwise isomerisation of the *trans*-isomer at higher temperature in a polar solvent.



The extent of isomerisation can be determined by various analytical techniques, for example by visible spectroscopy. *Cis*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride exists in form of two optical isomers (Λ -*cis* and Δ -*cis*) and their separation needs more demanding techniques.

Work

According to the following procedure prepare *cis*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride.

Chemicals

- *trans*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride, freshly prepared,
- ethanol, CH₃CH₂OH, denatured spirit,
- diethyl ether, (C₂H₅)₂O.

Procedure

In a small beaker dissolve 2,0 g of *cis*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride in 5 cm³ of water. Let stand the solution for 5 – 10 minutes. Put the beaker into a pre-heated water bath (Fig. 2). The solution changes its colour to violet. Evaporate *cca* half of its volume. Then, cool the mixture with ice. After cooling, filter out the violet *cis*-[CoCl₂(en)₂]Cl on a filter funnel with glass frit. Wash the product with small amount of cold ethanol and diethyl ether.




Fig. 2 Apparatus for preparation of *cis*-bis(ethane-1,2-diamine)-dichlorido-cobalt(III) chloride.

A note on isomerism: There are two structural (constitutional, geometrical) isomers (*cis*- and *trans*-) of the [CoCl₂(en)₂]⁺ cation. Additionally, two optical isomers (enantiomers, Λ -*cis*- and Δ -*cis*-) of the *cis*-isomer are known. If we write a chemical equation in which the *cis*-[CoCl₂(en)₂]Cl is a product, we always mean the mixture of two optical isomers, produced in statistically absolutely equal amounts. The final mixture containing the both enantiomers is called *racemic mixture* or *racemate*. A sample with only a single enantiomer is an *enantiomerically pure* or *enantiopure* compound.





Fig. 3 Visible spectra of isomerization *trans*-[CoCl₂(en)₂]Cl → *cis*-[CoCl₂(en)₂]Cl. Measured in methanol at 50 °C at the Graz University (Austria).

Safety instructions

Cobalt(II) chloride hexahydrate – $\text{CoCl}_2 \cdot 6 \text{H}_2\text{O}$

- R22** Harmful if swallowed.
R49 May cause cancer by inhalation.
R42/43 May cause sensitization by inhalation and skin contact.
R50/53 Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
S22 Do not breathe dust
S45 In case of accident or if you feel unwell seek medical advice immediately (show the label where possible).
S53 Avoid exposure – Obtain special instructions before use.
S60 This material and its container must be disposed of as hazardous waste.
S61 Avoid release to the environment. Refer to special instructions/safety data sheet.

Ethane-1,2-diamine (Ethylenediamine) – $\text{NH}_2\text{CH}_2\text{CH}_2\text{NH}_2$

- R10** Flammable.
R34 Causes burns.
R21/22 Harmful in contact with skin and if swallowed.
R42/43 May cause sensitization by inhalation and skin contact.
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S45 In case of accident or if you feel unwell seek medical advice immediately (show the label where possible).
S1/2 Keep locked up and out of the reach of children.
S36/37/39 Wear suitable protective clothing, gloves and eye/face protection.

Hydrochloric acid – HCl

- R34** Causes burns.
R37 Irritating to respiratory system.
S2 Keep out of the reach of children
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

Oxygen – O_2

- S21** When using do not smoke.

Dichlorido-bis(ethane-1,2-diamine)cobalt(III) chloride, $[\text{CoCl}_2(\text{en})_2]\text{Cl}$

- R36/37/38** Irritating to eyes, respiratory system and skin.
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S36 Wear suitable protective clothing.

Ethyl alcohol – $\text{C}_2\text{H}_5\text{OH}$

- R11** Highly flammable.
S7 Keep container tightly closed.
S16 Keep away from sources of ignition – No smoking.

Diethyl ether – $(\text{C}_2\text{H}_5)_2\text{O}$

- R12** Extremely flammable.
R19 May form explosive peroxides.
S9 Keep container in a well-ventilated place.
S16 Keep away from sources of ignition – No smoking.
S29 Do not empty into drains.
S33 Take precautionary measures against static discharges.