# Preparation of sodium aqua-(ethylenediaminetetraacetato)ferrate(III) dihydrate

Fig. 1 Ethylenediaminetetraacetic acid H<sub>4</sub>edta and its disodium salt Na<sub>2</sub>H<sub>2</sub>edta.

Ethylenediaminetetraacetate anion edta<sup>4-</sup> is an important hexadentate ligand. It forms very stable chelate complexes with many metal cations. The edta<sup>4-</sup> anion is coordinated to a central atom by four oxygen atoms and two nitrogen atoms (Fig. 2).

**Fig. 2** Structure of the complex anion Fe(edta)(H<sub>2</sub>O)]<sup>-</sup>.

Sodium aqua-(ethylenediaminetetraacetato)ferrate(III) dihydrate  $Na[Fe(edta)(H_2O)] \cdot 2H_2O$  may be prepared by two-step synthesis from the disodium salt  $Na_2H_2edta$ .

$$\begin{aligned} &H_2edta^{2-}(aq)+2\ OH^{-}(aq) \longrightarrow edta^{4-}(aq)+2\ H_2O(l) \\ &[Fe(H_2O)_6]^{3+}(aq)+edta^{4-}(aq) \longrightarrow [Fe(edta)(H_2O)]^{-}(aq)+5\ H_2O(l) \end{aligned}$$

#### Work

According to the procedure, prepare sodium aqua-(ethylenediaminetetraacetato)ferrate(III) dihydrate.

## **Chemicals**

- sodium dihydrogen ethylenediaminetetraacetate dihydrate Na<sub>2</sub>H<sub>2</sub>edta·2H<sub>2</sub>O, white crystalline substance,
- iron(III) chloride hexahydrate FeCl<sub>3</sub>·6H<sub>2</sub>O, dark yellow hygroscopic crystalline substance,
- sodium hydroxide, NaOH, white granules or flakes.

#### **Procedure**

In a 100 cm<sup>3</sup> beaker dissolve 0,4 g of NaOH in 10 cm<sup>3</sup> of water. Add 3,8 g of Na<sub>2</sub>H<sub>2</sub>edta·2H<sub>2</sub>O. Heat the mixture carefully until the solid substance dissolves. In a smaller beaker dissolve 2,5 g of FeCl<sub>3</sub>·6H<sub>2</sub>O in 5 cm<sup>3</sup> of water. Under continuous stirring, pour the prepared solution to the Na<sub>2</sub>H<sub>2</sub>edta and NaOH solution. Evaporate water by moderate boiling until the most of substance is excluded from the solution. Then, cool the reaction mixture in an ice bath and filter on the Büchner funnel. Wash the product thoroughly with ice water to remove all iron(III) cations. Finally, wash the product twice with ethanol and dry it by filter papers at room temperature.

### Safety instructions

## Sodium dihydrogen ethylenediaminetetraacetate dihydrate - Na<sub>2</sub>H<sub>2</sub>edta·2H<sub>2</sub>O

**R22** Harmful if swallowed.

R36/37/38 Irritating to eyes, respiratory system and skin.

**S2** Keep out of the reach of children.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**S36** Wear suitable protective clothing.

**S46** If swallowed, seek medical advice immediately and show this container or label.

## Iron(III) chloride hexahydrate - FeCl<sub>3</sub>·6H<sub>2</sub>O

R22 Harmful if swallowed.R38 Irritating to skin.

**R41** Risk of serious damage to eyes.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**S39** Wear eye/face protection.

### Sodium hydroxide – NaOH

R34 Causes burns.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S28 After contact with skin, wash immediately with plenty of (to be specified by the manufacturer).
S45 In case of accident or if you feel unwell seek medical advice immediately (show the label where

possible).

**S36/37/39** Wear suitable protective clothing, gloves and eye/face protection.